

Assay of Hydrogen Peroxide




H₂O₂: Mwt.=34.02 gm/mole

*Different H₂O₂ solution: 30%, 27%w/w and
6%, 3%w/v*


- H₂O₂ is slightly more viscous than water.
- Aqueous solution of H₂O₂ colorless and odorless. Rapidly decomposes when contact with organic matter and certain metals; Cu, Fe, Mn and if allowed to stay in alkaline solution.



- In acidic solution Fe^{2+} is oxidized to Fe^{3+} (H_2O_2 acting as an oxidizing agent):



- hydrogen peroxide acts as a reducing agent in its reaction with KMnO_4 *in our experiment*.




- **Uses:**

It is a *disinfectant* or *antiseptic* for treating wounds for its antimicrobial activity.


it is an effective *cleaning agent* prevents infection transmission in the hospital environment.

Procedure:

$\approx 0.1\text{N KMnO}_4$ 

10ml of unknown
25ml of D.W.
5ml of 50%v/v H₂SO₄ 

The buret contains the titrant.



This flask contains the solution to be titrated and the indicator.

End point: colorless  pink

- **Notes:**

1. We make the media acidic **to**;

- Prevent decomposition of H_2O_2

- Prevent reaction of $\text{KMnO}_4 \longrightarrow \text{MnO}_2$

- Increase oxidizing power of KMnO_4

2. Titration should be done slowly


- **Chemical principle:**



** You must write the oxidation and reduction equations to calculate the eqwt. for both KMnO_4 & H_2O_2 .*

- **Chemical factor:**

According to the chemical equation;



then complete the derivation of the Ch. Factor.

The Chemical factor is equal to:



- **Calculation:**
- Correct the **volume** of KMnO_4
 $N \times V = N' \times V' \text{ (corrected)}$
- Multiply the corrected volume by the Ch. Factor.
to get the **weight (wt.)** of H_2O_2
- the **%w/v** of $\text{H}_2\text{O}_2 = \text{wt.} / 10 \times 100$

- There is 2 methods to calculate the concentration of hydrogen peroxide in commercial products:

1. %W/V:

10ml of concentrated sol. was diluted to 200ml with D.W. \longrightarrow Then 20ml is taken from the dil. one and titrated with approximately 0.1N KMnO_4

<u>conc.</u>	<u>diluted</u>] $x = 20 \times 10 / 200 = \underline{1\text{ml}}$ of conc. sol. (original sol.)
10	200	
x	20	

Correct the volume of $\text{KMnO}_4 = \text{Vol.}$

$\text{Vol.} \times \text{ch. Factor} = \text{wt. of } \text{H}_2\text{O}_2 \text{ in } \underline{1\text{ml.}}$

Then calculate %w/v


2. Volume strength of the solution:

Decomposition occurs according to this equation:



Example: if the %w/v of H_2O_2 is 8.5%w/v *i.e.* It contains 8.5 gm of H_2O_2 in 100ml solution.

<u>wt.</u>	<u>Vol.</u>	
8.5	100ml	} X = 0.085 gm of H_2O_2 in 1ml
X	1ml (from the definition)	



So the volume strength of the sample is 28.