

**Preparation & Standardization of ≈ 0.1 N
Potassium Permanganate (KMnO_4) solution**

Redox Reaction



◦ *Preparation of $\approx 0.1N$ $KMnO_4$ solution:*

- Dissolve 3.3 gm of solid $KMnO_4$ in a 500 ml beaker containing water, stir thoroughly.
- Heat the solution on a hotplate for 15min. at $70-80^\circ C$.
- Cool the solution to room temp. and allow to stand for 2 days in a dark place.
- Filter the solution through the glass wool or asbestos into a 1L volumetric flask, complete the volume to the mark with D.W..



◦ **Standardization of KMnO_4 solution:**

- *Filter the prepared soln. through glass wool if it was left for a period of time (a week or more).*
- *Rinse and Fill the burette with filtered KMnO_4 soln.*
- *Weigh accurately about 200mg sod. oxalate ($\text{Na}_2\text{C}_2\text{O}_4$), previously dried at 110°C to constant weight.*
- *Dissolve it with 100ml of D.W. in a conical flask. Add 7ml. of conc. H_2SO_4 .*
- *Heat the soln. to about 70°C . (adjust temp. using a thermometer).*
- *Titrate against KMnO_4 soln. until a permanent pale pink color appears (persist for 15 sec.).*
- *Record the volume of KMnO_4 consumed.*



- **OXIDATION- REDUCTION REACTIONS:**

- ❖ Methods of analysis involve a change in valence of the reacting substances, thus there is always a transfer of electrons.
- ❖ Electrons lost by reducing agent are gained by oxidizing agent. The oxidizing agent will be reduced; gain of electrons and then decreasing in oxidation number, but the reducing agent will be oxidized; losing of electrons and then increasing in oxidation number.



- Chemical principle:
- Standardization of KMnO_4 against sod. oxalate follows oxidation-reduction reaction.
- Permanganate ion (MnO_4^-) acts as an oxidizing agent in acidic media, being reduced to Mn^{++} :



- Oxalate ion ($\text{C}_2\text{O}_4^{=}$) acts as reducing agent, being oxidized to CO_2 :



- KMnO_4 is widely used in volumetric analysis as a standard and oxidizing agent because:
 1. It is strong oxidizing agent.
 2. The intense color of KMnO_4 solution is sufficient to signal the end point in most titrations thus it is used as self indicator.
- On the other hand,
 - 1) KMnO_4 is a strong oxidizing agent, then it will oxidize any organic matter present in D.W., so heating for 15min. is to accelerate the oxidation of these organic matter.




- 2) KMnO_4 soln. is allowed to stand for 2 days to ensure the completion of the decomposition reaction.
- 3) It is filtered through asbestos to remove all traces of manganese dioxide (MnO_2) .
- 4) KMnO_4 is unstable in the presence of direct sunlight, organic matter ex. Filter paper.
- 5) H_2SO_4 is added since KMnO_4 is better oxidizing agent in acidic media, H_2SO_4 is added to keep $[\text{H}^+]$ ion conc. constant through out the titration.
- 6) Heat to about 70°C because the oxidation of sod. oxalate is rapid enough if the temp. above 60°C .



- Sod. oxalate, not oxalic acid, is used because:
 1. The salt can be obtained in very pure condition.
 2. Free from water of crystallization.
 3. It can be dried to 130 °c without decomposition.
 4. Stable to air.

The chemical equation for the std. of pot. Permanganate:







Calculations

$$\frac{\text{KMnO}_4}{N \times V} = \frac{\text{Na}_2\text{C}_2\text{O}_4}{\frac{\text{Wt.}}{\text{eq.wt.}} \times 1000}$$

67gm/mole

