University of Baghdad
College of Science for Women
Department of Chemistry



BIOCHEMISTRY LAB

(For Biology students/First class)

Prepared and Design by

Dr. Noor Ulhuda Ghazi Mohammed

Supervised by


Assist Professor Dr. Ahmed Younus

Assist Professor Dr. Israa Fadhil Ascar

2023-2024

5. Iodine Test: (General test of polysaccharides)

The iodine test is a common method for detecting the presence of polysaccharides, particularly starch, in a sample. The test is based on the interaction (adsorb I₂) between iodine and the helical structure of polysaccharides, which results in a color change when the iodine forms a complex with the polysaccharide.



Control Starch

-ve test +ve test

Method:


1 mL of the fresh starch solution & 1 drop of the iodine solution.

27/16 If the yellow color of the iodine reagent simply becomes diluted, no starch is present.

Record the observation as positive (blue) or negative (yellow).

Note: On heating (poly saccharide-Iodine) complexes losses blue colour because of dissociating the complex, while on cooling the blue colour again appears due to re-association of complex. Also Iodine test affected by base and acid medium. Below equation show when mineral acid is added to the basic medium, iodine is released again and a blue color appears.

$$3 I_2 + 6 NaOH \rightarrow 5 NaI + NaIO_3 + 3 H_2O$$



27/**18**