


Preparation and Standardization of 1 N NaOH Solution

- ❑ Sodium hydroxide is a strong base that is usually used to prepare standard alkaline solutions useful for volumetric analysis of acidic compounds.
- ❑ Sodium hydroxide is hygroscopic and can react with atmospheric carbon dioxide.



contaminant
(water soluble)



Preparation and Standardization of 1 N NaOH Solution

2. preparation of 1L of 1N NaOH solution

Dissolve 45 g sodium hydroxide in distilled water using a beaker, allow to cool, transfer the solution to a 1L volumetric flask and complete the volume to the mark by distilled water.

Then standardize against standard HCl solution.



Preparation and Standardization of 1 N NaOH Solution

□ standardization

The principle of the reaction is acid base titration reaction



- 1 N HCl solution is used as a secondary standard
- phenolphthalein is used as the indicator



Preparation and Standardization of 1 N NaOH Solution

□ procedure

- wash the burette with the D. W. and the titrant (NaOH)
- fill the burette with NaOH to a level (adjust it)
- wash a 10 mL – bulb pipette with D. W. then by a little of HCl solution; fill it to the mark with the acid
- transfer the acid into a clean conical flask; add D.W. (50 mL)
- add 2 drops of phenolphthalein indicator
- start titration by adding NaOH solution drop wise with continuous stirring until the solution changes from colourless to pink
- record the volume of NaOH solution used and calculate the normality


Note

wash the burette with water thoroughly




Preparation and Standardization of 1 N NaOH Solution

□ titration apparatus



Preparation and Standardization of 1 N HCl Solution

burette adjustment




Preparation and Standardization of 1 N HCl Solution

□ burette adjustment



How to read the burette?



Preparation and Standardization of 1 *N* NaOH Solution




**end point
(pink)**



Preparation and Standardization of 1 N NaOH Solution

□ calculations



N_1 :the normality of NaOH solution

V_1 :the volume of NaOH solution used

N_2 :the normality of HCl

V_2 : volume of HCl solution used (10 mL in our experiment)



Preparation and Standardization of 1 *N* NaOH Solution

□ Home work

Why have you used 45 g of NaOH to prepare 1000 mL of 1 *N* NaOH solution?

