# PREPARATION & STANDARDIZATION OF 0.1 N NAOH SOLUTION

Assist. Lec. Marwa Hashim

## **OBJECTIVE**

- - Prepare 0.1N NaOH solution
- - Standardize it using standard HCI

### PREPARATION PROCEDURE

- - Weigh 4g of NaOH pellets
- Dissolve in distilled water
- - Transfer to a 1L volumetric flask and fill to the mark

### STANDARDIZATION METHOD

- 1. Pipette 10mL of standard HCl into a conical flask
- 2. Add 1-2 drops of phenolphthalein indicator
- 3. Fill burette with prepared NaOH
- 4. Titrate until a faint pink color appears

# FINAL CALCULATION

• Use formula:  $N_1 \times V_1 = N_2 \times V_2$